

# Validation of cytology media for whole genome sequencing

<b>Functional Area</b>	<b>Cancer</b>	<b>Document Key</b>	PAR-SOP-004
<b>Document Owner</b>	Clare Craig	<b>Status</b>	FINAL
<b>Document Author</b>	Clare Craig	<b>Version</b>	1.0
<b>Document Reviewer(s)</b>	Shirley Henderson	<b>Version Date</b>	
	Kay Lawson	<b>Next Review Date</b>	
	Afshan Siddiq		
	Helen Stevens		
	Greg Elgar		
<b>Document Approval</b>	Tom Fowler		
	Sandi Deans		
<b>Electronic Signature</b>	<i>Tom Fowler</i> <a href="#">Tom Fowler (Feb 23, 2018)</a>	<b>Approval Date</b>	Feb 23, 2018
	<i>Sandi Deans</i> <a href="#">Sandi Deans (Feb 23, 2018)</a>		Feb 23, 2018
<b>Impact on Competent Personnel</b> (please choose Y/N in the boxes to the right)		<b>Read and understand</b>	Y/N
		<b>Re-train</b>	Y/N
<b>Transaction ID</b>	Adobe Sign Transaction Number: CBJCHBCAABAAPMy3ANGxleRN_UAix4_16sV5FaWAs2Uo		

Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

This document is uncontrolled if printed

## Contents

1	Document History and Control.....	3
1.1	Version History.....	3
2	Purpose & Scope.....	3
3	Roles & Responsibilities.....	3
4	Background to validation cohort.....	4
5	General Information.....	4
5.1	Consenting.....	5
5.2	Germline DNA sampling.....	5
5.3	Tumour sampling.....	5
5.4	Data entry.....	5
5.5	Sample submission.....	5
5.6	Return of Results.....	6
5.7	Where pilot cases differ from normal practice.....	6
6	Procedure.....	8
6.1	Process Flow.....	8
6.2	Procedural Steps.....	8
7	Additional Data.....	9
8	Definitions & Abbreviations.....	10
9	Data Model Requirements.....	10
10	DNA Extraction Protocol.....	10
11	Permission to Commence Collection of Cytology samples for the 100,000 genomes cancer programme.....	11
12	Notification.....	11
13	Related Documents, References & Procedures.....	11
14	Requirements for success.....	11
15	Appendix:.....	12

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 1 Document History and Control

The controlled copy of this document is maintained in the Genomics England internal document management system. Any copies of this document held outside of that system, in whatever format (for example, paper, email attachment), are considered to have passed out of control and should be checked for currency and validity. This document is uncontrolled when printed.

### 1.1 Version History

Version	Date	Description
0.1	04/12/2017	First draft
0.2	18/12/2017	Sandi and Shirley' edits
0.3	15/01/2018	Adding paragraph on permission to sample
0.4	29/01/2018	Incorporating Sandi Dean's edits
0.5	05/02/18	New title and checklist

## 2 Purpose & Scope

To validate the use of diagnostic fine needle aspiration (FNA) and fluid cytology samples for whole genome sequencing, this protocol has been devised.

NHS GMCs can opt to take part in a pilot of this technique and the requirements for this are outlined in this document.

## 3 Roles & Responsibilities

<i>Role</i>	<i>Responsibility</i>
Clinical Lead for validation cohort	Nominated NHS GMC individual responsible for ensuring the protocol is adhered to; samples are tracked and data is collected and submitted correctly.
Cytologist	The NHS GMC cytologist is responsible for diagnosis and tumour content assessment to identify eligible samples with sufficient residual material.
Scientist	The NHS GMC Scientist is responsible for keeping material fresh and judging when there is sufficient to enable the splitting of samples before preparation of cell blocks to retain fresh material for genomics.

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 4 Background to validation cohort

For many cancer patients a cytology sample will be the only diagnostic material available e.g. for metastatic effusions; pancreatic cancer; fine needle aspiration of metastatic carcinoma. Preliminary experimental work has indicated that DNA of sufficient quantity and quality for WGS is achievable with cytology samples. We would therefore like to initiate a pilot cohort at NHS GMCs to validate cytological samples for whole genome sequencing

This material is used for morphological assessment and diagnosis which often includes immunocytochemical staining. For the latter, the entire sample will be fixed in formalin which is detrimental to successful whole genome sequencing. Where there is sufficient material available consideration should be given where possible to dividing the sample and only fixing part of it for immunocytochemical staining and enabling fresh material to be kept separately for whole genome sequencing.

DNA yield from cytology samples must be sufficient for PCR free Whole Genome Sequencing. Samples with only enough DNA for PCR based library preparation cannot be submitted as part of the validation cohort.

Validation of this sampling method will require approximately 50 cases. The number required may increase during the course of the pilot if initial data reveals unexpected findings that require further work to fully understand. If there are inadequate samples to validate a particular cell preservation medium or a particular extraction technique then the validation cohort may be extended to cover more samples of that particular type.

It is hoped that validating cytology samples for Whole Genome Sequencing (WGS) will increase access to genomic testing for patients at diagnosis or where cytological material is readily available. Samples will be submitted as part of the main programme with a Whole Genome Analysis result returned to the GMC and these samples will count towards the recruitment trajectory and be eligible for funding.

## 5 General Information

When a sample is submitted as part of this validation cohort then the way the tumour is sampled will be affected as well as how the data is collected and submitted. There is a requirement to track these cases through the Genomics England service desk as outlined below.

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

**This document is uncontrolled if printed**

For all other aspects the usual guidance applies to how patients should be approached for consenting and how sampling should be carried out. The current sample handling guidance can be found here: <https://www.genomicsengland.co.uk/information-for-gmc-staff/sample-handling-guidance/> and guidance on consenting can be found as part of the Genomics Education Programme here: <https://www.genomicseducation.hee.nhs.uk/courses/courses/consent-ethics/>

### 5.1 Consenting

This work is considered service development work. The samples are therefore diagnostic samples. The decision on how best to sample a tumour will take place after arrival of the specimen in the pathology department. The decision on the best use of any DNA retrieved will be taken once the diagnosis and yield of DNA are known. Sampling for DNA can take place in the patient's best interests prior to specific consent for participation in the 100,000 Genomes Project in accordance with the Consensus Statement, which can be accessed here: <https://www.genomicsengland.co.uk/information-for-gmc-staff/cancer-programme/cancer-mdt-engagement-pack/>

### 5.2 Germline DNA sampling

A germline peripheral blood sample must be taken in an EDTA tube and the extracted germline DNA must be submitted with the tumour DNA sample to UKB for plating as set out in the current Sample Handling Guidance.

### 5.3 Tumour sampling

Only tumour that has been collected fresh can be submitted for WGS. The tumours must fulfil the eligibility criteria detailed here: <https://www.genomicsengland.co.uk/information-for-gmc-staff/cancer-programme/eligibility/>

### 5.4 Data entry

The patient must have a completed Registration file including consent information and a completed sample metadata file and Clinical Sample Test csv file (for tumour content assessment) prior to sample submission to UKB.

### 5.5 Sample submission

After DNA extraction, samples must be submitted to UKB following the recommendations of the current Sample Handling Guidance. Samples from this pilot should not be a clinical priority for fast or medium track slots, however, if they meet the volume and concentration requirements for fast track cases then they can be included on a fast track plate. Where possible samples submitted as part of the pilot should fulfil the fast track criteria in order to expedite the validation work.

#### Validation of cytology media for whole genome sequencing:


Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 5.6 Return of Results

Samples submitted as part of the validation cohort will be processed through the main programme and results will be returned to NHS GMCs via the usual route (cancer interpretation portal). As with all 100,000 Genomes Project cancer cases results should be interpreted and reported in accordance with the NHSE Cancer Validation and Reporting working group guidance document (in preparation).

## 5.7 Where pilot cases differ from normal practice


### Conventional process



### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## Validation pilot process




### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 6 Procedure

### 6.1 Process Flow



### 6.2 Procedural Steps

1. Where extra cytological material can be easily sampled and there is a high suspicion of malignancy, this should be sent for genomic sampling e.g. 100ml of malignant pleural effusion (note this is more than the routine 50ml).
2. All samples should be kept refrigerated to avoid DNA contamination and to inactivate nucleases.

#### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----



3. Where a diagnosis of cancer has been made or is suspected and a cell block is going to be prepared an assessment should be made as to whether there is adequate material to keep a proportion fresh for genomic sampling.
4. Once a diagnosis of cancer has been made and where there are >40% malignant nuclei in the sample then the cytopathologist should notify the Scientists that any residual material should be set aside for genomic testing.
5. If adequate fresh material is available the patient should be consented for inclusion in the 100,000 Genomes Project.
6. The genomic sample can be kept refrigerated for up to 72 hours prior to DNA extraction.

## 7 Additional Data

In addition to the submission of sample metadata the following data must be collected on each case within the pilot.

:

1. Participant ID
2. Disease subtype
3. Diagnostic morphology code
4. Collection media used, if any: PBS / RPMI / ThinPrep etc
5. Time refrigerated before freezing / extraction
6. Was the sample frozen?
7. Tumour content assessment
8. DNA yield
9. Local DNA QC metrics

Any NHS GMCs participating in the pilot must submit Standard Operating Procedures for the DNA extraction method they propose to use for cytology samples prior to commencing the pilot. This method must be used for all cases in the pilot. These methods must be submitted to the Clinical Lead for the validation cohort when requesting permission to participate in the pilot.

This data must be sent to the Clinical Lead for the diagnostic cytology sampling validation cohort at Genomics England, monthly or after every 10 cases have been sampled whichever is soonest. Genomics England will then assess the following outcomes for these samples:

1. DNA yield
2. Tumour purity
3. Coverage
4. Variant calling rate
5. Pick up rate for known mutations including smoking related genomic signatures

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 8 Definitions & Abbreviations

Abbreviation / Term	Description
FNA	Fine Needle Aspirate
DNA	Deoxyribonucleic Acid
PBS	Phosphate Buffered Saline
RPMI	Roswell Park Memorial Institute cell culture medium
QC	Quality control
UKB	United Kingdom Biobank
WGS	Whole Genome Sequencing
NHS GMC	National Health Service Genomic Medicine Centre
FF	Fresh Frozen

## 9 Data Model Requirements

There are a number of specific requirements around entering sample metadata on these samples since the current data model does not fully support cytology samples.

Clinic Sample Type should be entered as DNA FF Tumour since there is no alternative enumeration for cytology samples in the current data model. Laboratory Sample ID from the FluidX tube for the FNA DNA sample should be documented.

Tissue Source should be entered as 'non\_standard\_biopsy'.

## 10 DNA Extraction Protocol

All cytology samples and particularly fine needle aspiration samples can have cohesive sheets of cells so it is important that the homogenisation and lysate steps are completed as for any fresh tumour sample.

The DNA extraction technique used is critical to achieving adequate DNA yield and quality. Extraction must be performed using a kit and protocol which has been recommended by the manufacturer as appropriate for cells pelleted from cell suspension and using a suitable volume. QC requirements of the resulting DNA are the same as for a fresh frozen tissue sample.

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 11 Permission to Commence Collection of Cytology samples for the 100,000 genomes cancer programme

Prior to commencing collection of cytology samples, NHS GMCs will be required to have submitted a completed checklist (see appendix) and reviewed procedural aspects of the cytology sample collection and the details of the DNA extraction method they plan to use with Genomics England and NHS England (usually via a teleconference).

## 12 Notification

In order that we may carefully track and monitor the metadata and sequencing metrics on these samples the Genomics England Helpdesk ([ssd.genomics@hscic.gov.uk](mailto:ssd.genomics@hscic.gov.uk)) must be notified of each case with the subject title 'Cytology Sample', the DNA extraction technique used and both the Participant ID and Laboratory Sample ID (for the DNA samples) should be included in the text of the notification.

## 13 Related Documents, References & Procedures

Sampling of tumours with this technique should adhere to the guidance set out in the current [Sample Handling Guidance](#).

## 14 Requirements for success

Cytology samples will allow genomic sequencing from cancer patients who otherwise could not get a sequence. The standard that these samples need to reach in order to be of diagnostic benefit therefore needn't be as high as the gold standard. For example, if purity or coverage were not as good as for fresh frozen samples but the samples were still able to pick up mutational signatures with an acceptable degree of sensitivity then there will still be diagnostic value in sequencing these samples.

Results of the validation from each GMC cohort will be shared with the participating GMC once they are available.

Once results are available for the validation cohort as a whole these will be shared with all GMCs.

### Validation of cytology media for whole genome sequencing:

Document Key	PAR-SOP-004	Version Number	1.0
--------------	-------------	----------------	-----

## 15 Appendix:

Checklist for participation in cytology media validation pilot	
<b>Sample Handling</b>	
Proposed tumour types that will be sampled	
Number of cases expected	
Confirmation of which laboratories will process the samples for DNA extraction	
Storage and transportation requirements prior to DNA extraction	
Details of proposed DNA extraction method and SOP	
<b>Data capture</b>	
Confirmation of how Pathology Department will record cases including data on ineligible cases and those that fail QC <ul style="list-style-type: none"> <li>i) Tissue sampling details</li> <li>ii) DNA extraction details</li> </ul>	
Designated responsibility for collecting the meta-data	
<b>Sample tracking</b>	
Process for tracking when these patients have been diagnosed, consented to GEL and sample has been retrieved, including responsibility to complete this action	
Confirmed process to notify GEL service desk that a sample has been submitted and associated patient identifier, including responsibility to complete this action	
<b>Teleconferences</b>	
Proposed operators to talk through intended sampling methods	
Clinical Scientist to talk through DNA extraction methods proposed	